
Scientific Abstracts

10 years of clinical research | **Vol. 1**

Dear Doctor,

The present collection of scientific abstracts gives you an overview of the existing scientific evidence for mectron PIEZOSURGERY®. Ten years ago, Dr. Tomaso Vercellotti came up with the idea of piezoelectric bone surgery. From the beginning of this revolutionary surgical technique, many universities and specialized surgeons have been involved in the clinical and scientific validation of the piezoelectric bone surgery. Today, a continuously growing number of publications confirm the clinical benefits of mectron PIEZOSURGERY® :

- micrometric cuts
- selective cuts
- maximum intra-operative visibility

With the scientific proof statement of these benefits, mectron PIEZOSURGERY® has defined a new dimension in bone surgery and today is considered worldwide the only evidence based alternative to traditional surgical tools.

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Wallace SS, Mazor Z, Froum SJ, Cho SC, Tarnow DP.
 Schneiderian membrane perforation rate during sinus elevation using piezosurgery: clinical results of 100 consecutive cases.
 Int J Periodontics Restorative Dent. 2007; 27(5):413-419 10

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Vercellotti T, Russo C, Gianotti S
 A New Piezoelectric Ridge Expansion Technique in the Lower Arch - A Case Report.
 World Dentistry. 2000 online article 12

Enislidis G, Wittwer G, Ewers R.
 Preliminary Report on a Staged Ridge Splitting Technique for Implant Placement in the Mandible: A Technical Note.
 Int J Oral Maxillofacial Implants. 2006; 21(3):445-449 12

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	<p>Stübinger S, Robertson A, Zimmerer SK, Leiggenger C, Sader R, Kunz C. Piezoelectric Harvesting of an Autogenous Bone Graft from the Zygomaticomaxillary Region: Case Report. Int J Periodontics Restorative Dent. 2006; 26(5):453-457</p>	16
	<p>Sivolella S, Berengo M, Scarin M, Mella F, Martinelli F. Autogenous particulate bone collected with a piezo-electric surgical device and bone trap: a microbiological and histomorphometric study. Arch Oral Biol. 2006; 51(10):883-891</p>	16
	<p>Happe A. Use of a piezoelectric surgical device to harvest bone grafts from the mandibular ramus: report of 40 cases. Int J Periodontics Restorative Dent. 2007; 27(3):241-249</p>	17
	<p>Patel A, Schofield J. Using Piezosurgery to harvest a block bone from the symphyseal region: a clinical case presentation. Implant Dentistry Today. 2007; 1(4):20-24</p>	17
	<p>Sohn DS, Ahn MR, Lee WH, Yeo DS, Lim SY. Piezoelectric osteotomy for intraoral harvesting of bone blocks. Int J Periodontics Restorative Dent. 2007; 27(2):127-131</p>	19
	<p>Landes CA, Stübinger S, Laudemann K, Rieger J, Sader R. Bone harvesting at the anterior iliac crest using piezoosteotomy versus conventional open harvesting: a pilot study. Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2008; 105(3):e19-e28</p>	20
IMPLANT SITE PREPARATION	<p>Preti G, Martinasso G, Peirone B, Navone R, Manzella C, Muzio G, Russo C, Canuto RA, Schierano G. Cytokines and Growth Factors Involved in the Osseointegration of Oral Titanium Implants Positioned using Piezoelectric Bone Surgery Versus a Drill Technique: A Pilot Study in Minipigs. J Periodontol. 2007; 78(4):716-722</p>	21
OSTEOTOMY CLOSE TO NERVES	<p>Bovi M. Mobilization of the Inferior Alveolar Nerve with simultaneous implant insertion: A New Technique. A Case Report. Int J Periodontics Restorative Dent. 2005; 25(4):375-383</p>	22
	<p>Geha H, Gleizal A, Nimeskern N, Beziat JL. Sensitivity of the Inferior Lip and Chin following Mandibular Bilateral Sagittal Split Osteotomy Using Piezosurgery. Plast Reconstr Surg. 2006; 118(7):1598-1607</p>	23
	<p>Schaeren S, Jaquiéry C, Heberer M, Tolnay M, Vercellotti T, Martin I. Assessment of nerve damage using a novel ultrasonic device for bone cutting. Département de Parodontologie de l'UFR de Stomatologie et Chirurgie Maxillo-Faciale. J Oral Maxillofac Surg. 2008; 66(3):593-596</p>	24
	<p>Sakkas N, Otten JE, Gutwald R, Schmelzeisen R. Transposition of the mental nerve by piezosurgery followed by postoperative neurosensory control: a case report. Br J Oral Maxillofac Surg. 2008; 46(4):270-271</p>	25

EXTRACTIONS	Grenga V, M. Bovi. Piezoelectric Surgery for Exposure of Palatally Impacted Canines. J Clin Orthod. 2004; 38(8):446-448	25
	Sivolella S, Berengo M, Fiorot M, Mazzuchin M. Retrieval of blade implants with piezosurgery: two clinical cases. Minerva Stomatol. 2007; 56(1-2):53-61	26
PERIODONTAL SURGERY	Vercellotti T, Pollack AS. A New Bone Surgery Device: Sinus Grafting and Periodontal Surgery. Compend Contin Educ Dent. 2006; 27(5):319-325	27
DISTRACTION OSTEOGENESIS	González-García A, Diniz-Freitas M, Somoza-Martín M, García-García A. Piezoelectric Bone Surgery Applied in Alveolar Distraction Osteogenesis: A Technical Note. Int J Oral Maxillofac Implants. 2007; 22(6):1012-1016	27
	Menini I, Zornitta C, Menini G. Distraction Osteogenesis for Implant Site Development Using a Novel Orthodontic Device: A Case Report. Int J Periodontics Restorative Dent. 2008; 28(2):189-196	28
HISTOLOGICAL ADVANTAGES	Vercellotti T, Nevins ML, Kim DM, Nevins M, Wada K, Schenk RK, Fiorellini JP. Osseous Response following Resective Therapy with a Piezosurgery®. Int J Periodontics Restorative Dent. 2005; 25(6):543-549	28
	Gleizal A, Li S, Pialat JB, Béziat JL. Transcriptional expression of calvarial bone after treatment with low - intensity ultrasound: An in vitro study. Ultrasound Med Biol. 2006; 32(10):1569-1574	29
TECHNIQUE ADAVANTAGES	Vercellotti T. La Chirurgia Ossea Piezoelettrica. Il Dentista Moderno. 2003; 5:21-55	30
	Boioli LT, Vercellotti T, Tecucianu JF. La chirurgie piézoélectrique: Une alternative aux techniques classiques de chirurgie osseuse. Inf Dent. 2004; 86(41):2887-2893	30
	Vercellotti T. Technological characteristics and clinical indications of piezoelectric bone surgery. Minerva Stomatol. 2004; 53(5):207-21	32
	Lambrecht JT. Intraorale Piezo-Chirurgie. Schweiz Monatsschr Zahnmed. 2004; 114(1):29-36	33
	Siervo S, Ruggli-Milic S, Radici M, Siervo P, Jäger K. Piezoelektrische Chirurgie. Eine alternative Methode für die schonende Chirurgie Schweiz Monatsschr Zahnmed. 2004; 114(4):365-377	34

Troiani C, Russo C, Ballarani G, Vercellotti T. Piezoelectric Surgery: A new reality to cut and manage bone. Maxillo Odontostomatologia- International Journal of Maxillo Odontostomatology - S.I.M.O. 2005; 4(1):23-28	34
Vercellotti T. Características tecnológicas e indicaciones clinicas de la cirugia òsea piezoelétrica. Revista Mundo Dental. 2005; 26-28	35
Schlee M. Ultraschallgestützte Chirurgie - Grundlagen und Möglichkeiten. Z Zahnärztl Impl. 2005; 21(1):48-59	36
Schlee M, Steigmann M, Bratu E, Garg AK. Piezosurgery: basics and possibilities. Implant Dent. 2006; 15(4):334-340	37
Sembronio S, Albiero AM, Polini F, Robiony M, Politi M. Intraoral endoscopically assisted treatment of temporomandibular joint ankylosis: Preliminary report. Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2007; 104(1):e7-e10	37
Su YC. Development and clinical application of ultrasonic osteotomy in dentistry. Shanghai kou qiang yi xue = Shanghai journal of stomatology. 2007; 16(1):1-7 [chinese]	38
Peivandi A, Bugnet R, Debize E, Gleizal A, Dohan DM. Piezoelectric osteotomy: applications in periodontal and implant surgery. Rev Stomatol Chir Maxillofac. 2007; 108(5):431-440	39
Stübinger S, Landes C, Seitz O, Zeilhofer HF, Sader R. Ultraschallbasiertes Knochenschneiden in der Oralchirurgie: eine Übersicht anhand von 60 Patientenfällen. Ultraschall Med. 2008; 29(1):66-71	39
Labanca M, Azzola F, Vinci R, Rodella LF Piezoelectric surgery: Twenty years of use. Br J Oral Maxillofac Surg. 2008; 46(4):265-269	41
Stübinger S, Filippi A, Sader R, Zeilhofer HF. Intraoral Piezosurgery: Preliminary Results of a New Technique. J Oral Maxillofac Surg. 2005; 63(9):1283-1287	41

Vercellotti T, De Paoli S, Nevins M.

The Piezoelectric Bony Window Osteotomy and Sinus Membrane Elevation: Introduction of a New Technique for Simplification of the Sinus Augmentation Procedure.

ENGLISH

Int J Periodontics Restorative Dent 2000; 20(4):359-365

■ The Piezoelectric Bony Window Osteotomy and Sinus Membrane Elevation: Introduction of a New Technique for Simplification of the Sinus Augmentation Procedure.

All of the surgical techniques to elevate the maxillary sinus present the possibility of perforating the schneiderian membrane. This complication can occur during the osteotomy, which is performed with burs, or during the elevation of the membrane using manual elevators. The purpose of this article is to present a new surgical technique that radically simplifies maxillary sinus surgery, thus avoiding perforating the membrane. The piezoelectric bony window osteotomy easily cuts mineralized tissue without damaging the soft tissue, and the piezoelectric sinus membrane elevation separates the schneiderian membrane without causing perforations. The elevation of the membrane from the sinus floor is performed using both piezoelectric elevators and the force of a physiologic solution subjected to piezoelectric cavitation. Twenty-one piezoelectric bony window osteotomy and piezoelectric sinus membrane elevations were performed on 15 patients using the appropriate surgical device (Mectron Piezosurgery System). Only one perforation occurred during the osteotomy at the site of an underwood septa, resulting in a 95% success rate. The average length of the window was 14 mm; its height was 6 mm, and its thickness was 1.4 mm. The average time necessary for the piezoelectric bony window osteotomy was approximately 3 minutes, while the piezoelectric sinus membrane elevation required approximately 5 minutes.

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Piezoelectric Bone Surgery for Sinus Bone Grafting.

ENGLISH

The Sinus Bone Graft, Second Edition. Edited by Ole T. Jensen, Quintessence Books. 2006; 23:273-279

■ Piezoelectric Bone Surgery for Sinus Bone Grafting.

Over the last half century, applications of ultrasonic instrumentation in bone surgery largely failed, particularly in light of reports of their substantial ineffectiveness or even negative influence on the healing process. Studies comparing different types of instruments (ultrasonic instrument, bur, and

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scalpel) have obtained conflicting results based on a variety of experimental protocols.

Several authors explored the possibility of using ultrasonics in maxillary sinus surgery, but the effort was curtailed because of insurmountable difficulties. Significant limiting factors include the failure of traditional ultrasonic instruments to cut bone structures that are more than 1 mm thick and a tendency for excessive overheating leading to necrosis. However, recent scientific-technologic research of bone's response to mechanical stimuli, particularly ultrasonic vibration, has led to a new and highly effective technique for osteotomy and osteoplasty called piezoelectric bone surgery.

The revolutionary properties of piezoelectric surgery have simplified many common osseous surgical procedures, including sinus bone grafting.

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■ Contribution of Piezosurgery to sinus elevation

Summary: Piezoelectric surgery is a new technique providing an interesting alternative to traditional bone surgery techniques. The insert performs a selective cutting action on mineralised tissues but is ineffective on the soft tissues. It can be used for osteotomies and for osteoplastic surgery exploiting the microvibrations of the insert. Because of the small size of the insert, it is possible to obtain cuts of sub-millimetric width. Furthermore, the cooling of the insert by the physiological liquid contributes towards achieving a bloodless surgical field. The advantages of piezosurgery by comparison with the techniques currently used for sinus elevation are presented below.

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Wallace SS, Mazor Z, Froum SJ, Cho SC, Tarnow DP.

Schneiderian membrane perforation rate during sinus elevation using piezosurgery: clinical results of 100 consecutive cases.

ENGLISH

Int J Periodontics Restorative Dent. 2007; 27(5):413-419

■ Schneiderian membrane perforation rate during sinus elevation using piezosurgery: clinical results of 100 consecutive cases.

The lateral window sinus elevation procedure has become a routine and highly successful preprosthodontic procedure that is used to increase bone volume in the posterior maxilla for the placement of dental implants. Many surgical techniques have been proposed that provide access to the maxillary sinus through the lateral wall to allow for elevation of the sinus membrane. Among these are the multiple variations of the hinge and complete osteotomy techniques, which make use of rotary cutting instruments for the antrostomy. The most common intraoperative complication with these surgical approaches is perforation of the schneiderian membrane, with perforation rates of 14% to 56% reported in the literature. In most instances, perforation occurs either while using rotary instruments to make the window or when using hand instruments to gain initial access to begin the elevation of the membrane from the sinus walls. This article presents an alternative approach that uses a piezoelectric instrument for the sinus elevation procedure. Although new to the United States, this approach has been used successfully in Europe for many years. The membrane perforation rate in this series of 100 consecutive cases using the piezoelectric technique has been reduced from the average reported rate of 30% with rotary instrumentation to 7%. Furthermore, all perforations with the piezoelectric technique occurred during the hand instrumentation phase and not with the piezoelectric inserts.

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■ Piezoelectric Surgery. Maxillary sinus lift techniques.

Bone surgery, the one used for the maxillary sinus lift in Implantology in particular, has to follow two main purposes. The first purpose consists in the predictability of the results which are wished to be obtained by means of surgery so that a suitable therapeutic valence could be reached, according to the clinical needs of each patient. The second

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purpose consists in reaching the goals just described before with the less morbidity as possible during the post-operative course, both in the sinus region where bone grafting is performed and in the potential donor site from which the graft is harvested (Chanavaz 1990; Sinus Graft Consensus Conference 1996). The purpose of this chapter is to introduce the new Piezoelectric Surgery techniques dedicated to sinus lift pointing out the surgeon intraoperative safety and the patient lowered morbidity.

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Vercellotti T.

Piezoelectric Surgery in Implantology: A Case Report - A New Piezoelectric Ridge Expansion Technique.

ENGLISH Int J Periodontics Restorative Dent. 2000; 20(4):359-365

■ Piezoelectric Surgery in Implantology: A Case Report - A New Piezoelectric Ridge Expansion Technique.

The purpose of this preliminary article is to present a new surgical technique that, thanks to the use of modulated-frequency piezoelectric energy scalpels, permits the expansion of the ridge and the placement of implants in single-stage surgery in positions that were not previously possible with any other method. The technique involves the separation of the vestibular osseous flap from the palatal flap and the immediate positioning of the implant between the 2 cortical walls. The case report illustrates the ridge expansion and positioning of implants step by step in bone of quality 1 to 2 with only 2 to 3 mm of thickness that is maintained for its entire height. To obtain rapid healing, the expansion space that was created for the positioning of the implant was filled, following the concepts of tissue engineering, with bioactive glass synthetic bone graft material as an osteoconductive factor and autogenous platelet-rich plasma as an osteoinductive factor. The site was covered with a platelet-rich plasma membrane. A careful evaluation of the site when reopened after 3 months revealed that the ridge was mineralized and stabilized at a thickness of 5 mm and the implants were osseointegrated.

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Vercellotti T, Russo C, Gianotti S.

A New Piezoelectric Ridge Expansion Technique in the Lower Arch - A Case Report.

ENGLISH

World Dentistry. 2000; online article

■ A New Piezoelectric Ridge Expansion Technique in the Lower Arch - A Case Report.

This article presents a new piezoelectric ridge expansion technique that permits the use of the bone flap procedure even in areas of highly mineralized bone, such as in the anterior region of the mandible, thus avoiding the risk of fracture. The new piezoelectric ridge expansion technique makes use of a specifically designed piezoelectric device with bone scalpels that vibrate at a variable ultrasound frequency, cutting bone easily, precisely, and non-traumatically.

The case reported is that of a patient with an edentulous mandibular ridge defect with thickness reduced to 3 mm. The bone flap expansion technique was performed, followed by a piezoelectric ridge osteotomy to cut the ridge and obtain separation of the buccal and palatal osseous flaps, allowing the positioning of the implants in a single stage surgery. The defect obtained by the distraction of the two cortical walls was filled with bioactive glass mixed with platelet concentrate (PC) and covered with a platelet rich plasma (PRP) membrane. When the site was reopened after six months, the ridge bone appeared to be well stabilized with a thickness of 5 mm, an increase of 2 mm; the implants also appeared to be osteointegrated.

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Enislidis G, Wittwer G, Ewers R.

Preliminary Report on a Staged Ridge Splitting Technique for Implant Placement in the Mandible: A Technical Note.

ENGLISH

Int J Oral Maxillofac Impl. 2006; 21(3):445-449

■ Preliminary Report on a Staged Ridge Splitting Technique for Implant Placement in the Mandible: A Technical Note.

Purpose: Narrow edentulous alveolar ridges less than 5 mm wide require horizontal augmentation for the placement of screw-type dental implants. A staged approach to ridge splitting in the mandible to decrease the risk of malfracture during osteotomy is presented.

Materials and Methods: Five consecutive patients with 6 long-span edentulous areas of the mandibular ridge were included in this study. After

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corticotomy of a rectangular buccal segment and a 40-day healing period, the mandibular ridge was split, leaving the buccal periosteum attached to the lateralized segment. Seventeen dental implants were placed, and the gap between the implants and the bone filled with a mixture of venous blood and a porous algae-derived hydroxyapatite.

Results: All buccal segments fractured as planned at the basal corticotomy during ridge splitting. After 6 months, all implants were stable and surrounded by bone; prosthetic loading with fixed partial dentures was successful in all cases.

Discussion: In the mandible, greenstick fracture during widening with osteotomes has not been controllable to date because of cortical thickness of the bone; the risk of malfracture during single-stage ridge splitting was high. With this approach, the location of the greenstick fracture is predetermined, and the perfusion for the buccal segment remains intact, although vascularization shifts from internal perfusion from spongy bone after the first intervention to external perfusion from the periosteum after the second intervention. The buccal cortical segment remains a pedicled graft after ridge splitting.

Conclusion: The preliminary results of this report indicate that staged ridge splitting can be a safe technique which overcomes the problems associated with single-stage ridge expansion/ridge splitting procedures without causing significant delay in treatment.

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Vercellotti T, Podestà A.

Orthodontic Microsurgery: A New Surgically Guided Technique for Dental Movement.

ENGLISH

Int J Periodontics Restorative Dent. 2007; 27:325-331

■ Orthodontic Microsurgery: A New Surgically Guided Technique for Dental Movement.

Eight patients with malocclusions were treated with a new orthodontic-surgical technique that reduces the duration of treatment compared to conventional techniques. The monocortical tooth dislocation and ligament distraction (MTDLD) technique combines two different dental movements that work separately but simultaneously on opposite root surfaces. On the root surface corresponding to the direction of movement, vertical and horizontal microsurgical corticotomies are performed around each tooth root with a piezosurgical microsaw to eliminate cortical bone resistance. The immediate applica-

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tion of strong biomechanical forces produces rapid dislocation of the root and the cortical bone together. On the root surface opposite the direction of movement, the force of dislocation produces rapid distraction of ligament fibers. During the osteogenic process that follows, application of normal orthodontic biomechanics achieves the final tooth movement. All eight patients underwent periodontal and radiologic examinations for more than 1 year after treatment. No periodontal defects were observed in any of the patients, including one with a severe malocclusion and a thin periodontal tissue biotype. Compared to traditional orthodontic therapy, the average treatment time with the MTDLD technique in the mandible and maxilla was reduced by 60% and 70%, respectively.

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Chiriac G, Herten M, Schwarz F, Rothamel D, Becker J.

Autogenous bone chips: influence of a new piezoelectric device (Piezosurgery®) on chips morphology, cell viability and differentiation.

ENGLISH

J Clin Periodontol. 2005; 32(9):994-999

■ Autogenous bone chips: influence of a new piezoelectric device (Piezosurgery®) on chips morphology, cell viability and differentiation.

Aim: The aim of the present study was to investigate the influence of a new piezoelectric device, designed for harvesting autogenous bone chips from intra-oral sites, on chip morphology, cell viability and differentiation.

Methods: A total of 69 samples of cortical bone chips were randomly gained by either (1) a piezoelectric device (PS), or (2) conventional rotating drills (RD). Shape and size of the bone chips were compared by means of morphometrical analysis. Outgrowing osteoblasts were identified by means of alkaline phosphatase activity (AP), immunohistochemical staining for osteocalcin (OC) synthesis and reverse transcriptase-polymerase chain reaction phenotyping.

Results: In 88.9% of the RD and 87.9% of the PS specimens, an outgrowth of adherent cells nearby the bone chips was observed after 6-19 days. Confluence of cells was reached after 4 weeks. Positive staining for AP and OC identified the cells as osteoblasts. The morphometrical analysis revealed a statistically significant more voluminous size of the particles collected with PS than RD.

Conclusion: Within the limits of the present study, it may be concluded that both the harvesting methods are not different from each other concerning their

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detrimental effect on viability and differentiation of cells growing out of autogenous bone chips derived from intra-oral cortical sites.

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■ Piezosurgery and pre-implant preparation of bone - Grafting by adding blocks of autogenous bone with portions of harvested bone.

Abstract: Piezoelectric surgery is a new technique that is an interesting alternative to traditional bone surgery techniques. The insert performs a selective cutting action on mineralised tissues, while it is ineffective on the soft tissues. It can be used to carry out osteotomies and osteoplasties thank to the micro-vibrations of the active part of the inserts. Because of the small size of the insert, it is possible to make sub-millimetric cuts. Furthermore, the fact that the insert is cooled by physiological liquid contributes towards keeping the surgical field bloodless. Two clinical cases of rehabilitation of vestibulo-palatal bone defects are presented below. Two autogenous bone grafts, taken from the mandibular retromolar triangle, were carried out using the piezosurgery technique to reconstruct the bone and create a volume sufficient to receive an implant.

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Piezoelectric Harvesting of an Autogenous Bone Graft from the Zygomaticomaxillary Region: Case Report.

ENGLISH

Int J Periodontics Restorative Dent. 2006; 26(5):453-457

■ Piezoelectric Harvesting of an Autogenous Bone Graft from the Zygomaticomaxillary Region: Case Report.

The zygomaticomaxillary region offers a large amount of cortical bone that can be obtained simply and safely using the precise and selective cutting properties of a piezosurgical device. A block from this area fits nicely into anterior or premolar maxillary recipient sites and is thus the ideal choice, as no secondary surgical field is needed. As in conventional sinus bone graft procedures, the complication rate is minimal and after a 5-month healing period, the augmented region can be used for stable and esthetic oral implant placement.

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Sivolella S, Berengo M, Scarin M, Mella F, Martinelli F.

Autogenous particulate bone collected with a piezo-electric surgical device and bone trap: a microbiological and histomorphometric study.

ENGLISH

Arch Oral Biol. 2006; 51(10):883-891

■ Autogenous particulate bone collected with a piezo-electric surgical device and bone trap: a microbiological and histomorphometric study.

The aims of this study were to determine the microbiological and particle size characteristics of particulate bone collected with a piezosurgery device and bone trap, and to reduce bacterial contamination after treatment of debris with rifamycin SV. Samples were taken from 10 patients who underwent surgical extraction of their third lower molars. The ostectomy was performed with a piezosurgery device, and the debris was collected with a surgical aspiration set equipped with a bone trap. Two aliquots were taken from each sample, one of which was treated with rifamycin SV. The second aliquot, used as a control, was treated with a physiological solution. In the samples immersed in antibiotic solution, there was a statistically significant ($P < 0.005$) reduction in bacterial contamination. The stringent protocol followed in this study has proved valid for collection of material, and treatment with rifamycin SV was found to reduce bacterial contamination in collected material.

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Happe A.

Use of a piezoelectric surgical device to harvest bone grafts from the mandibular ramus: report of 40 cases.

ENGLISH

Int J Periodontics Restorative Dent. 2007; 27(3):241-249

■ Use of a piezoelectric surgical device to harvest bone grafts from the mandibular ramus: report of 40 cases.

Forty patients (12 men, 28 women; mean age 44.9 +/- 13.95 years) with inadequate bone volume for implant placement in either the maxilla or the mandible were scheduled for bone-grafting procedures. Bone grafts were harvested with a piezoelectric surgical device (Piezosurgery, Mectron) from 45 donor sites and grafted into recipient sites prior to implant placement. After healing of the grafted sites, 109 implants were placed. Clinical evaluations were performed to assess the size and quality of the obtained grafts, complications at the donor and recipient sites, integration of the grafts, bone quality, resorption of the grafts, and ability to properly place implants. Piezosurgery allowed precise, clean, and smooth cutting with excellent visibility. Mean graft size was 1.15 cm³ (SD 0.5), with a maximum of 2.4 cm³. The quality of the obtained grafts was mainly cortical (71%), 42 of the 45 donor sites healed uneventfully (93.33%), and 50 of the 52 grafted sites healed uneventfully (96.15%).

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Patel A, Schofield J.

Using Piezosurgery to harvest a block bone from the symphyseal region: a clinical case presentation.

ENGLISH

Implant Dentistry Today. 2007; 1(4):20-24

■ Using Piezosurgery to harvest a block bone from the symphyseal region: a clinical case presentation.

The placement of titanium implants in sites that are deficient of alveolar ridge width has always been unpredictable due the lack of bone around the implant. Implants placed into the available bone can often lead to a compromised aesthetic result due to poor positioning and direction of the implants, especially in the aesthetic zone. A number of surgical techniques have been proposed over the years to create adequate bone volume. These include the use of xenograft, alloplast and allograft materials with the use of resorbable or non-resorbable membranes (Garg, 1999). It has been shown that autogenous bone grafts are still the most predictable method to achieve good bone quality and quantity for the

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implant site (Garg, 1999).

Donor sites for autogenous bone include extraoral and intraoral sites. Extraoral sites allow for harvesting of large amounts of bone. These sites involve complex surgical procedures that usually require general anaesthesia with an increase in morbidity at the donor sites (Chiapasco et al, 1999, Sindet-Pedersen and Enemark, 1990). It has been shown that in the facial skeleton, membranous bone (such as that grafted from the mandible) undergoes less resorption than endochondral bone, such as the iliac crest (Zins and Whitacker, 1979).

The main advantage of an intraoral donor site is that harvesting of bone can be performed under local anaesthesia. There are two main intraoral donor sites for autologous bone grafting available to the implantologist (Triplett and Silhow, 1998). The choice of site is based on the size of the defect to be augmented and the type and quantity of bone needed and the access to the donor site (Raghoobar et al, 2001). Only a certain amount of bone can be harvested from the retromolar area. It has been observed that the amount of bone volume obtained from the retromolar area is half of what can be harvested from the symphyseal region (Misch 1997). Retromolar grafts are usually used for single tooth augmentations.

Following harvesting of symphysis bone grafts, some long-term neurosensory disturbances have been described (Nkene et al, 2001, Raghoobar et al, 2001). In a prospective study evaluating intraoral donor sites for bone grafting, it was reported that 29% of symphysis graft patients had altered sensitivity of their mandibular incisor teeth (Misch 1997). All findings were resolved within six months. The majority of studies illustrate a continuous improvement of tooth sensitivity over time. The authors have audited 94 consecutive symphysis grafts performed over the last 10 years. All patients have reported some altered sensation in their lower incisors. The altered sensation lasted between three and six months in all but two cases. Of these two cases, one patient had altered sensation for 11 months. The other patient, on whom a very large chin graft was harvested, had some discoloration of her lower incisors. Despite this, sensitivity testing revealed normal results. The discoloration disappeared after four months. In contrast, surgical access to the retromolar area might be compromised because of anatomy, i.e. the location of the mandibular canal and width of the ramus. This donor site appears to be associated with less permanent neurosensory disturbances (Nkenke et al, 2002).

When using conventional bone cutting instruments such as reciprocating microsaws or trephines, there is a risk of postoperative necrosis. Comparative histological studies have shown necrosis on cut bone caused by surgical burs, trephines and microsaws (Aro et al, 1981). Over the past five years, piezo-

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surgery has been introduced as a new technique for osteotomy and osteoplasty using ultrasonic vibration. The piezosurgery blades allow for maximum intraoperative precision and minimal tissue damage, and the selective frequency of the scalpel minimises the risk to adjacent soft tissues. Histological studies have also shown minimal cellular damage to the resected bone margins, bony matrix and underlying marrow spaces. A recent study concluded that piezosurgery provided more favourable osseous repair and remodelling than with conventional surgical burs with surgical osteotomy and osteoplasty procedures (Vercellotti et al, 2005).

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Sohn DS, Ahn MR, Lee WH, Yeo DS, Lim SY.

Piezoelectric osteotomy for intraoral harvesting of bone blocks.

ENGLISH

Int J Periodontics Restorative Dent. 2007; 27(2):127-131

■ Piezoelectric osteotomy for intraoral harvesting of bone blocks.

Grafting with intraoral bone blocks is a good way to reconstruct severe horizontal and vertical bone resorption in future implant sites. The Piezosurgery System (Mectron) creates an effective osteotomy with minimal or no trauma to soft tissue, in contrast to conventional surgical burs or saws. In addition, piezoelectric surgery produces less vibration and noise because it uses microvibration, in contrast to the macrovibration and extreme noise that occur with a surgical saw or bur. Microvibration and reduced noise minimize a patient's psychologic stress and fear during osteotomy under local anesthesia. The purpose of this article is to describe the harvesting of intraoral bone blocks using the piezoelectric surgery device.

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Bone harvesting at the anterior iliac crest using piezoosteotomy versus conventional open harvesting: a pilot study.

ENGLISH

Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2008; 105(3):e19-e28

■ Bone harvesting at the anterior iliac crest using piezoosteotomy versus conventional open harvesting: a pilot study.

Objective: The objective of this study was to evaluate piezoosteotomy for bone harvesting at the iliac crest in a pilot collective.

Study design: Thirteen patients underwent anterior iliac crest bone graft harvesting using piezoosteotomy. These patients were compared to an age- and sex-matched retrospective cohort of 13 patients who underwent bone graft harvesting at the anterior iliac crest using conventional chisels and saws in an otherwise identical protocol.

Results: Harvested bone volumes and operation times were comparable; hospitalization time was briefer and pain levels at the first 2 postoperative days less in the piezoosteotomy group; pain medication requirement was comparable. Spearman's-rho correlation showed a strong tie between pain medication requirement and harvested bone volume in the piezoosteotomy group and between harvested bone volume and operation time in the conventional surgery collective.

Conclusions: The correlation between required pain medication and harvested bone volume together with lower pain levels in the piezoosteotomy group indicate pain more exclusively related to local osseous damage than to correlated soft tissue lacerations. Conventional technique shows a correlation between operation time and harvested bone volume indicating a more time-effective procedure, although total operation time was comparable. Piezoosteotomy appears to have potential use in iliac crest bone harvesting.

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Cytokines and Growth Factors Involved in the Osseointegration of Oral Titanium Implants Positioned using Piezoelectric Bone Surgery Versus a Drill Technique: A Pilot Study in Minipigs.

ENGLISH

J Periodontol. 2007; 78(4):716-722

■ **Cytokines and Growth Factors Involved in the Osseointegration of Oral Titanium Implants Positioned using Piezoelectric Bone Surgery Versus a Drill Technique: A Pilot Study in Minipigs.**

Background: Most dental implants are positioned using a drilling surgery technique. However, dentistry recently experienced the implementation of piezoelectric surgery. This technique was introduced to overcome some of the limitations involving rotating instruments in bone surgery. This study used biomolecular and histologic analyses to compare the osseointegration of porous implants positioned using traditional drills versus the piezoelectric bone surgery technique.

Methods: Porous titanium implants were inserted into minipig tibias. Histomorphology and levels of bone morphogenetic protein (BMP)-4, transforming growth factor (TGF)-beta2, tumor necrosis factor alpha, and interleukin-1beta and -10 were evaluated in the peri-implant osseous samples.

Results: Histomorphological analyses demonstrated that more inflammatory cells were present in samples from drilled sites. Also, neo-osteogenesis was consistently more active in bone samples from the implant sites that were prepared using piezoelectric bone surgery. Moreover, bone around the implants treated with the piezoelectric bone surgery technique showed an earlier increase in BMP-4 and TGF-beta2 proteins as well as a reduction in proinflammatory cytokines.

Conclusion: Piezoelectric bone surgery appears to be more efficient in the first phases of bone healing; it induced an earlier increase in BMPs, controlled the inflammatory process better, and stimulated bone remodeling as early as 56 days post-treatment.

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Bovi M.

Mobilization of the Inferior Alveolar Nerve with simultaneous implant insertion: A New Technique. A Case Report.

ENGLISH

Int J Periodontics Restorative Dent. 2005; 25(4):375-383

■ Mobilization of the Inferior Alveolar Nerve with simultaneous implant insertion: A New Technique. A Case Report.

For reconstruction of an atrophied posterior mandible, different therapeutic options have been proposed, such as autologous bone grafting, guided bone regeneration for vertical ridge augmentation, and inferior alveolar nerve (IAN) mobilization with simultaneous implant placement. The possible dehiscence of soft tissues covering the surgical zone makes the first and second techniques unpredictable. Moreover, two surgical sites are necessary and a long treatment time is required (about 12 months). With IAN mobilization, only one surgical intervention is required and the total treatment time is shorter (about 6 months). However, this technique risks irreversible damage to the IAN, with consequent functional alterations. Current studies have shown extreme variability in the examination of functionality of the neurovascular bundle after its mobilization. This variability can be attributed both to the methodology used for the tests, which evoke subjective answers from the patient, and to the surgical procedure, which is highly dependent on operator technique. Nerve damage can be the result of an overstretched mucoperiosteal flap in the premolar area to achieve optimal visibility of the surgical zone. This article reports a case in which a new surgical approach to IAN mobilization is performed using a specifically engineered device for simplified bone surgery (Mectron Piezosurgery). This device enables the surgeon to cut hard tissue without injuring the soft tissues. Therefore, there is a lower risk of damaging the IAN, and it is possible to reduce overstretching of the mental nerve by creating a smaller bone window and using an apico-coronal inclination of instruments to capture the neurovascular bundle.

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Geha H, Gleizal A, Nimeskern N, Beziat JL.

Sensitivity of the Inferior Lip and Chin following Mandibular Bilateral Sagittal Split Osteotomy Using Piezosurgery.

ENGLISH

Plast Reconstr Surg. 2006; 118(7):1598-1607

■ Sensitivity of the Inferior Lip and Chin following Mandibular Bilateral Sagittal Split Osteotomy Using Piezosurgery.

Background: Bimaxillary osteotomy, including bilateral sagittal split osteotomy, is the most commonly performed orthognathic surgical intervention in Europe and the United States. Neurosensory perturbation in the territory of the inferior alveolar nerve is a reported adverse effect of bilateral sagittal split osteotomy. Piezosurgery is a relatively new technique that allows bone to be cut while preserving soft tissues, including nerves. The purpose of this study was to assess inferior alveolar nerve function through clinical neurosensory testing after bilateral sagittal split osteotomy using Piezosurgery.

Methods: Between February and September of 2004, 20 patients (40 sides) presenting with dento-skeletal deformities underwent bimaxillary osteotomy, including bilateral sagittal split osteotomy. The Mectron Piezosurgery device was used to perform all sagittal splits, with distraction being performed between the two bone valves. The inferior alveolar nerve was evaluated both objectively with clinical neurosensory testing, including pin-prick sensation, light touch sensation, and two-point discrimination tests, and subjectively at the following time points: preoperatively; at days 5, 7, and 10, postoperatively; and at the second month postoperatively.

Results: The anatomical integrity of the inferior alveolar nerve was respected in all cases. Observed normal results for the different tests at 10 days were 90, 82, and 70 percent, respectively, for pin-prick sensation, light touch sensation, and two-point discrimination. By computing scores for different clinical neurosensory tests, the authors observed between 75 and 80 percent complete neurosensory recuperation as early as the second postoperative month.

Conclusions: Piezosurgery used for bilateral sagittal split osteotomy allows prompt recovery of inferior alveolar nerve neurosensory function within 2 months. No comparison is possible with the results using the standard technique for bilateral sagittal split osteotomy.

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Schaeren S, Jaquiéry C, Heberer M, Tolnay M, Vercellotti T, Martin I.
Assessment of nerve damage using
a novel ultrasonic device for bone cutting.

ENGLISH

J Oral Maxillofac Surg. 2008; 66(3):593-596

■ **Assessment of nerve damage using a novel ultrasonic device for bone cutting.**

In oral, maxillofacial, and spinal surgery, osteotomy often must be performed in close vicinity to nerve tissue, with the potential risk of transient or permanent neurologic injuries (eg, to trigeminal nerve branches). Traditional tools, such as rotating burs and oscillating saws, are highly effective in cutting bone tissue but are not selective for bone, and thus can produce significant harm to surrounding soft tissues, especially nerves.

Piezosurgery (Mebiotec, Sestri Levante, Italy), a device for bone-cutting based on low-frequency (25 to 29 kHz) ultrasonic vibrations, has recently been introduced in oral and maxillofacial surgery. This device can improve the control and precision of osteotomy and improve bone healing due to reduced local trauma. A recent in vitro study of the use of Piezosurgery for transposition of the inferior alveolar nerve in the mandibles of cadaver sheep showed that this technique caused roughening of the epineurium without affecting deeper structures and induced less injury than a conventional rotary bur. A pilot clinical study found that Piezosurgery reduced neurosensory disturbances in orthognathic surgery of the mandible.

To the best of our knowledge, no experimental in vivo study has yet demonstrated that bone cutting using Piezosurgery prevents damage to soft tissue, especially nerves. The aim of the present study was to assess the potential damage of Piezosurgery to a peripheral nerve on direct contact, in 2 possible scenarios.

The first scenario corresponds to the situation in which the surgeon does not immediately realize the contact with the nerve and continues to apply the same force used for cutting bone (working force) for an estimated reaction time of 5 seconds. The second scenario corresponds to an accidental slip of the device on the bone tissue, causing a direct hit on the nerve with a higher force (peak force), estimated to be twice the working force, but for a shorter time (1 second). The rat sciatic nerve was selected for the study because it is an established model for studying sensory and motor nerve damage and recovery and also is of equivalent size of branches of human trigeminal nerve.

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Sakkas N, Otten JE, Gutwald R, Schmelzeisen R.

Transposition of the mental nerve by piezosurgery followed by postoperative neurosensory control: a case report.

ENGLISH

Br J Oral Maxillofac Surg. 2008; 46(4):270-271

■ Transposition of the mental nerve by piezosurgery followed by postoperative neurosensory control: a case report.

Transposition of the mental nerve is a preprosthetic procedure that is effective for patients with hyperaesthesia caused by the effect of a dental prosthesis on the alveolar ridge. We present the case of a 74-year-old woman with pain and hyperaesthesia of the right inferior alveolar nerve caused by a dental prosthesis. Caudal transposition of the right mental nerve by piezosurgery resulted in postoperative neurosensory controls of the lower lip showing normal nerve function 2 months later.

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■ ■ ■ EXTRACTIONS

Grenga V, M. Bovi.

Piezoelectric Surgery for Exposure of Palatally Impacted Canines.

ENGLISH

J Clin Orthod. 2004; 38(8):446-448

■ Piezoelectric Surgery for Exposure of Palatally Impacted Canines.

A precise surgical cut is extremely important when the crown of a palatally impacted maxillary canine is close to the roots or crowns of the central and lateral incisors, as often occurs. Traditional burs do not distinguish between the mineralization and hardness of the bone, the radicular cementum and the enamel. Piezoelectric surgery is precise enough to account for such differences, thus avoiding damage to the adjacent teeth.

The EX1 insert for the Mectron Piezosurgery hand-piece provides a good demarcation of the exposed crown. It also removes the follicle faster than manual instruments can, with no loss of cervical bone or consequent gingival recession.

The piezoelectric instrument controls bleeding during the surgical procedure, ensuring a dry field for bonding to the impacted tooth and eliminating the need for gauze or electrocoagulation. Total surgical time is greatly reduced with this method.

The piezoelectric device is quieter than conventional

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handpieces, and it avoids the trauma associated with a manual hammer and chisel. Therefore, it improves both the comfort and the cooperation of the patient.

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Sivolella S, Berengo M, Fiorot M, Mazzuchin M.

Retrieval of blade implants with piezosurgery: two clinical cases.

ENGLISH

Minerva Stomatol. 2007; 56(1-2):53-61

■ Retrieval of blade implants with piezosurgery: two clinical cases.

In this work an ultrasound device was used to perform an ostectomy for the removal of blade implants in order to save as much bone tissue as possible, so that root form implants might later be inserted. Two patients underwent surgery for the removal of two blade implants (one maxillary, the other mandibular) that were no longer functional. The peri-implant ostectomy was carried out with a piezoelectric surgery device. The instrument demonstrated to be effective and precise during ostectomy, providing an extremely thin cutting line. During the course of the operation and at controls after 7 and 30 days, patients did not show any relevant complications and both still had sufficient alveolar bone to be treated with root form implants. The piezosurgery device proved to be an effective instrument in interventions requiring a significant saving of bone tissue, extreme precision in cutting, and respect of soft tissues

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Vercellotti T, Pollack AS.

A New Bone Surgery Device: Sinus Grafting and Periodontal Surgery.

ENGLISH

Compend Contin Educ Dent. 2006; 27(5):319-325

■ A New Bone Surgery Device: Sinus Grafting and Periodontal Surgery.

The PiezoSurgery device is a new instrument developed specifically for osseous surgery that has applications in various dental and medical surgical specialties. Using low-frequency piezoelectric ultrasonic vibration, the PiezoSurgery device precisely cuts bone without cutting soft tissue. This article presents 3 cases illustrating the use of the device in sinus grafting and periodontal osseous surgery.

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Piezoelectric bone surgery is based on ultrasonic vibration of a device functioning as an osteotome. This device allows precise cuts to be made in bone structures without provoking lesions of adjacent soft tissues and at the same time offering excellent visibility within the surgical field. The use of this technique in alveolar distraction osteogenesis is described. Piezoelectric surgery appears to make the cutting of the transport segment easier and safer.

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Menini I, Zornitta C, Menini G.

Distraction Osteogenesis for Implant Site Development Using a Novel Orthodontic Device: A Case Report.

ENGLISH

Int J Periodontics Restorative Dent. 2008; 28(2):189-196

■ Distraction Osteogenesis for Implant Site Development Using a Novel Orthodontic Device: A Case Report.

The case report describes the distraction osteogenesis of a small bone fragment supporting a maxillary central incisor that had become ankylosed and infraoccluded after trauma. After initial orthodontic treatment, surgery was performed, entailing a trapezoidal osteotomy with minimal periosteal elevation through a horizontal incision. The osteotomy was performed using a piezoelectric surgical device (Piezosurgery, Mectron) with the thinnest inserts. On this occasion, a device was applied, anchored to the tooth, and activated by a micrometer screw. On completion of distraction, the fragment was retained for 5 months. The ankylosed tooth was then removed, and an implant was inserted and provided with an immediate provisional prosthesis without occlusal loading. Five months later, a definitive prosthesis with occlusal loading was applied and is still in place 1 year later without any problems.

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■■■ HISTOLOGICAL ADVANTAGES

Vercellotti T, Nevins ML, Kim DM, Nevins M, Wada K, Schenk RK, Fiorellini JP. Osseous Response following Resective Therapy with a Piezosurgery®.

ENGLISH

Int J Periodontics Restorative Dent. 2005; 25(6):543-549

■ Osseous Response following Resective Therapy with a Piezosurgery®.

A piezoelectric instrument vibrating in the ultrasonic frequency range was investigated for its potential use in periodontal resective therapy. The rate of postoperative wound healing (baseline and 14, 28, and 56 days after surgery) in a dog model following surgical ostectomy and osteoplasty was the marker used to compare the efficacy of this instrument (PS) with a commonly used carbide bur (CB) or a diamond bur (DB). The surgical sites treated by CB or DB lost bone, in comparison to baseline measurements, by the 14th day, while the surgical sites

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treated by PS revealed a gain in the bone level. By day 28, the surgical sites treated by all three instruments demonstrated an increased bone level and regeneration of cementum and periodontal ligament. However, by day 56, the surgical sites treated by CB or DB evidenced a loss of bone, versus a bone gain in the PS-treated sites. Thus, it appears that PS provided more favourable osseous repair and remodeling than CB or DB when surgical ostectomy and osteoplasty procedures were performed. Therefore, PS could be regarded as being efficacious for use in osseous surgery.

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Gleizal A, Li S, Pialat JB, Béziat JL.

Transcriptional expression of calvarial bone after treatment with low - intensity ultrasound: An in vitro study.

ENGLISH

Ultrasound Med Biol. 2006; 32(10):1569-1574

■ Transcriptional expression of calvarial bone after treatment with low - intensity ultrasound: An in vitro study.

The capacity to reossify a calvarial bone defect is very small in mature animals and in humans greater than 2 y of age. The clinical treatment of injured tissue sites of bones by low-intensity pulsed ultrasound is widespread, but little is known about the precise effects of ultrasound on the fundamental processes that promote repair and regeneration. In this study, we used real-time polymerase chain reaction (RT-PCR) to investigate the expression of osteogenesis-associated genes after stimulation by low-intensity ultrasound in adult mouse osteoblast from the parietal calvaria. The gene associated with the Runx2 pathway had notably higher levels after 1, 2 and 3 days of stimulation. Therefore, low-intensity ultrasound seems to have an effect on the transcriptional gene expression of the calvarial bone in vitro.

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Il Dentista Moderno. 2003; 5:21-55

■ Piezoelectric bone surgery.

This article aims at illustrating the origin, development and clinical advantages of piezoelectric surgery, a new, revolutionary bone cutting technique. First, the features of traditional instruments are described highlighting their limitations in terms of intraoperative precision and safety. Piezoelectric surgery, however, goes beyond these limits thanks to the physical characteristics of the mechanical action leading to osteotomy. Piezoelectric cutting can be defined as micrometric and selective in nature for it is not visible to the human eye (the cutting mechanical action is obtained thanks to the modulated, ultrasonic micro-vibration of only 60-200 mm of amplitude); it is selective for piezoelectric ultrasounds used at low frequency allow to cut mineralized tissue with great effectiveness, whilst they are inert on soft tissue. Thanks to these and other physical features, piezoelectric cutting offers a number of clinical advantages, such as cutting precision, safe performance and excellent intraoperative visibility. A number of histological results are also presented from the wealth of multidisciplinary scientific research that has characterized the development and evolution of this young surgical technique from the start. Dental indications are presented, such as, for example, piezoelectric impacted or ankylosed root extraction techniques, and socket debridement for immediate post-extraction implant placement. Lastly, operative protocols for periodontal surgery, ridge expansion and maxillary sinus lift are also detailed.

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■ **Piezoelectric surgery. An alternative to conventional bone-surgery techniques.**

Piezoelectric bone surgery, known as piezosurgery, is a technique that can be used to carry out osteotomies and osteoplasties thanks to the ultrasonic micro-vibrations of the inserts. This technique has been devised by Dr Vercellotti in order to broaden the indications of bone surgery, to simplify surgical protocols and improve the predictability of the results.

At the moment, bone burs and saws are used for this type of surgery. The burs are operated by micromotors and the user must counter the reactive movements induced by the torque of the instrument, thus forcing him to make an additional effort. This lowers the surgeon's sensitivity, above all when structures with different degrees of mineralization are encountered.

Similarly, saws produce macrovibrations that require resistance to be exerted while the cut is being made. Although saws ensure perfect linearity, they do not enable control of the depth of the cut, and this means that it is necessary to use additional instruments such as hammers and chisels to complete the incision.

When using the piezosurgery technique, on the other hand, the effort required to make a cut is very slight. This means that greater precision is achieved, guaranteed by the microvibrations of the insert. The thickness of the cut obtained will be determined solely by the thickness of the insert used, and the surgeon's macromovements will be minimal. The frequencies produce an effective cutting action only on the mineralised tissues, which makes for safer bone surgery, even in the proximity of anatomical structures to be preserved, such as vascular and nerve masses or other soft tissues. It is

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thus possible to minimise the risk of lesions to these structures, considering that they are not sensitive to the frequencies used. This gives rise to a high-precision technique, with a selective cutting action.

A histological study has enabled the healing mechanisms of the tissues following piezoelectric cutting to be studied. This study, carried out on dogs, showed that it was possible to note the absence of necrosis traces on the surface of the cut, and the presence of nucleated osteocytes, indicative of the only slightly traumatic effect of this technique. Before using piezoelectric surgery it is essential to acquire the technique and to become familiar with the instruments. The surgeon's movements differ from those typical of other bone surgery techniques, and the cut is generated by the microvibrations of the insert. It is necessary to find and acquire a perfect balance between the pressure exerted and the speed of movement of the insert.

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Vercellotti T.

Technological characteristics and clinical indications of piezoelectric bone surgery.

ENGLISH

Minerva Stomatol. 2004; 53(5):207-214

■ Technological characteristics and clinical indications of piezoelectric bone surgery.

Piezoelectric bone surgery, also simply known as piezosurgery, is a new technique for osteotomy and osteoplasty utilizing an innovative ultrasonic surgical apparatus. This technique was created and developed in response to the need to reach major levels of precision and safety in bone surgery, as compared to that available by the usual manual and motorized instruments. The instrument which offers these results, known as the Mectron Piezosurgery Device, is characterized by piezoelectric ultrasonic vibrations of a frequency of 29 kHz and a range between 60/200 Hz. Due to their characteristics, the microvibrations allow a selective cut of only mineralized structures without damaging soft tissues, which remain undamaged even in case of accidental contact. The micrometric vibration ensures precise cutting action and at the same time maintains a blood-free site because of the physical phenomenon of cavitation. The micrometric vibration makes the instrument manageable and permits major interoperative control with a consequent increase in safety especially in anatomically difficult areas. The absence of macrovibrations makes the instrument more manageable and allows greater interoperative control with a significant increase in the cutting safety in the more difficult anatomical cutting zones. Given its innovative nature, piezosurgery distinguishes itself as being different from the conventional techniques utilised in bone

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surgery. Consequently, it requires different surgical skill. To master the right surgical skill it is essential to demonstrate an adequate learning curve.

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■ Intraoral Piezo-Surgery

The Piezo „effect“ describes a physical interaction, which is used in different fields of clinical medicine. An ultrasound technology results in the power of a surgical cut technique, which is able to distinguish between hard and soft tissue. With the presented instruments it is possible to treat bone mechanically in a piezo-surgical manner without hurting the surrounding soft tissues when touching them.

Four examples:

- Tooth extraction before implantation under minimal grinding of bone,
- Gaining bone for peri-implant transplantation,
- Preparation of Sinus Lift,
- Exposing the inferior alveolar nerve

demonstrate the use of Piezo-Surgery in the intraoral area.

This method shows a remarkable addition to the intraoral operation techniques.

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■ Piezoelectric surgery. An alternative method of minimally invasive surgery

The piezoelectric surgery is an ultrasonic surgery which represents a novel and alternative method to the conventional hard and soft tissue management with rotating instruments in the oral district. The innovation of this technique is mainly due to three major characteristics: a) a micrometric cut with an amplitude between 60 µm and 200 µm, b) a selective cut which works on hard tissues but not on soft tissues, c) a relative blood free surgical field due to the air-water cavitation effect of the ultrasonic device. The final result of these characteristics is a definite clinical advantage with regard to the cut precision, the sparing of vital nervous soft tissues, the better visualisation of the surgical area. The present work shows through a series of case reports the possible use of the piezoelectric surgery in the oral district, presenting observed advantages and disadvantages. The entire clinical work was performed with the Piezosurgery device Mectron (Mectron Medical Technology, I-16042 Carasco).

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Troiani C, Russo C, Ballarani G, Vercellotti T.

Piezoelectric Surgery: A new reality to cut and manage bone.

ENGLISH

Maxillo Odontostomatologia- International Journal of Maxillo Odontostomatology - S.I.M.O. 2005; 4(1):23-28

■ Piezoelectric Surgery: A new reality to cut and manage bone.

Summary: Piezosurgery is a new osteotomic and osteoplastic technique that uses ultrasound technology as its functioning principle: this system aims to improve the precision and safety of the traditional instruments normally used in bone surgery. Its special characteristics include: the selective

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cutting due to the use of low frequency ultrasound of 29.000 Hz (respecting soft tissues, vessels, nerves, mucosa); the precision of the cut guaranteed by the micro vibrations of the insert (40-200 µm); the cavitation phenomenon which allows a blood free and clean site, absence of macro vibrations that guarantees greater control and greater operative safety even in the most complex anatomical zones. The surgeon with minor experience in different surgical techniques will need a reasonable amount of practice to achieve a satisfactory learning curve.

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■ Technological characteristics and clinical indications of the piezoelectric bone surgery.

Summary: Piezoelectric bone surgery, also known as Piezosurgery, is a new osteotomy and osteoplasty technique using an innovative piezoelectric surgery device. The technique in question was developed by the author in order to meet the need to achieve a higher level of precision and safety compared to manual and motorised instruments normally used in bone surgery. The main feature of the device which has enabled these results to be achieved, known as the Mectron Piezosurgery Device, is its piezoelectric vibration at a frequency of 29kHz and with an amplitude comprised between 60 and 200 µm. Thanks to their features, the microvibrations enable selective cutting of the mineralised structure without causing any damage to the soft tissues even if these are touched by accident. In addition, the

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micrometric vibrations enable great precision of cutting and at the same time they keep the surgical site clean (bloodless) thanks to the physical phenomenon of cavitation. The absence of macrovibrations makes it easy to grip the device in such a way as to enable greater intra-operative control, thus improving safety and precision of cutting. This means that it can also be used in anatomical areas that are difficult to access. Piezoelectric bone surgery is an innovative technique, and its use differs from traditional techniques. This means that the manual skills required are also different. They can be learnt according to a suitable learning curve.

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■ Piezosurgery - rudiments and possibilities

This article describes the therapeutical potential of the Piezoelectric Bone Surgery in implantology. PIEZOSURGERY enables more precise and non traumatic cutting of bone in comparison to conventional methods (micrometric cut). The instrument vibrates with a modulated ultrasonic frequency. Because the vibration frequency of Piezosurgery is optimal for mineralized tissue it does not cut soft tissue. In fact, in order to cut soft tissue other ultrasonic frequencies are necessary. This allows an osteotomy which prevents delicate anatomic soft tissue structures like nerves or the Schneiderian Membrane of being damaged (selective cut). The cooling liquid flows laminar along the tip. Thereby, a more effective cooling also in deeper regions is achieved. Methods like crest expansion, sinus lift elevation, harvesting of bone chips and bone blocks are becoming much more secure and easier.

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ENGLISH Implant Dent. 2006; 15(4):334-340

■ **Piezosurgery: basics and possibilities.**

Useful in a variety of oral surgery procedures, piezosurgery has therapeutic features that include a micrometric cut (precise and secure action to limit tissue damage, especially to osteocytes), a selective cut (affecting mineralized tissues, but not surrounding soft tissues), and a clear surgical site (the result of the cavitation effect created by an irrigation/cooling solution and oscillating tip). Because the instrument's tip vibrates at different ultrasonic frequencies, since hard and soft tissues are cut at different frequencies, a „selective cut“ enables the clinician to cut hard tissues while sparing fine anatomical structures (e.g., schneiderian membrane, nerve tissue). An oscillating tip drives the cooling-irrigation fluid, making it possible to obtain effective cooling as well as higher visibility (via cavitation effect) compared to conventional surgical instruments (rotating burs and oscillating saws), even in deep spaces. As a result, implantology surgical techniques such as bone harvesting (chips and blocks), crestal bone splitting, and sinus floor elevation can be performed with greater ease and safety.

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Sembronio S, Albiero AM, Polini F, Robiony M, Politi M.
Intraoral endoscopically assisted treatment of temporomandibular joint ankylosis: Preliminary report.

ENGLISH Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2007; 104(1):e7-e10

■ **Intraoral endoscopically assisted treatment of temporomandibular joint ankylosis: Preliminary report.**

Temporomandibular joint (TMJ) ankylosis is characterized by the formation of a bony or fibrous mass that replaces the normal articulation. To avoid a possible re-ankylosis it is mandatory to perform a radical, complete resection of the bony/fibrous mass. We treated a patient affected by right temporomandibular joint ankylosis performing the osteotomy of the ankylotic mass through a preauricular and intraoral approach under endoscopic control. The use of a piezoelectric instrument, whose microvibrations are adjusted to target only mineralized tissue, reduced the risk of damage to soft tissues. Then a temporalis muscle and fascia flap were used as the interpositional material. Through the endoscope it was easy to check the medial aspect of the resection

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and suture the flap. At 1-year follow-up the patient had significantly increased maximal mouth opening. No evidence of relapse of the joint ankylosis was shown by radiological studies. Intraoral endoscopic assistance may be useful to make the removal of the ankylotic mass safer, and the anchorage of the temporalis muscle and fascia flap more accurate, reducing the risk of re-ankylosis.

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Su YC.

Development and clinical application of ultrasonic osteotomy in dentistry.

ENGLISH

Shanghai kou qiang yi xue = Shanghai journal of stomatology. 2007; 16(1):1-7 [chinese]

■ Development and clinical application of ultrasonic osteotomy in dentistry.

Ultrasonic osteotomy is a new technique for osseous surgery utilizing an innovative ultrasonic surgical apparatus. Oral surgeons have started using this newly developed ultrasonic bone scalpel when performing precise osteotomies. This technique is known as piezosurgery also and has found wide applications in dental and medical specialities. Numerous protocols including conventional diamond round bur, safe scraper, rongeur pliers and gouge shaped bone chisel have been proposed to simplify osteotomy and osteoplasty in bone surgery before this technique was created. But only ultrasonic osteotomy could reach high levels of precision and safety as compared to those instruments. Ultrasonic osteotomy device is characterized by piezoelectric ultrasonic vibrations of a frequency of 25-29 kHz and a range between 60-200 µm. Due to low-frequency ultrasonic vibration, the device precisely cuts bone without cutting soft tissue even in case of accidental contact. It also maintains a blood-free site during cutting bone. This technique was first used in sinus floor elevation for bony window osteotomy and sinus membrane elevation. After several years, this technique is utilized in more and more dental indications such as ridge expansion, maxillary sinus floor elevation, harvesting of bone graft, extraction and socket debridement for the installation of immediate-loaded implants, extraction of third molars, periodontal surgery, cyst exeresis and inferior alveolar nerve transposition. It's also reported in hand surgery, cranial osteoplasty and segmental maxillary Le Fort osteotomy. Ultrasonic osteotomy is different from the conventional techniques utilized and proved to be a safe and useful technique in bone surgery.

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Peivandi A, Bugnet R, Debize E, Gleizal A, Dohan DM.

Piezoelectric osteotomy: applications in periodontal and implant surgery.

ENGLISH

Rev Stomatol Chir Maxillofac. 2007; 108(5):431-440

■ Piezoelectric osteotomy: applications in periodontal and implant surgery.

Reverse piezoelectricity is currently used in Dentistry for ultrasonic scaling devices. Using this technology, increasingly more powerful ultrasonic surgical lancets were developed within a few years, and these new tools have provided many practical solutions in oral and maxillofacial surgery. This article reviews the basic principles of piezosurgery, and its numerous clinical applications: in oral surgery (atraumatic extractions, wisdom teeth exposure, periapical cyst debridement, pre-prosthetic surgery), more specifically in periodontal (root planning and bone remodeling, crown lengthening) and implant surgery (sinus lift, removal of fractured implants, bone ridge augmentation, bone graft harvesting (crestal, retromolar, chin or tori)). The recent increased power of these instruments allows emphasizing a much wider range of clinical applications, extending to all surgical fields.

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■ **Ultrasonic Bone Cutting in Oral Surgery: a Review of 60 Cases.**

Purpose: Surgical reconstruction of bony defects in the oral cavity can often be challenging since thin and fragile bony structures are especially prone to fracture caused by bulky cutting tips or the application of significant pressure by conventional mechanical instruments. The risk of accidental damage to adjacent soft tissue structures, such as nerves, by a dental drill or saw is also extremely high. The use of modulated ultrasound (piezosurgery) makes it possible to overcome such complications as a result of the precise and minimally invasive surgery technique which is limited to mineralized hard tissue.

Materials and methods: In 60 patients (38 male, 22 female) a piezosurgery device was used for different bone augmentation procedures before dental implant placement. The instrument uses modulated ultrasound (25 - 30 kHz) and the amplitude of the working tip ranges from 60 µm to 200 µm. The device was employed for sinus floor elevation (25), alveolar ridge augmentation using an autogenous block graft (25), alveolar ridge splitting (5) or lateralization of the alveolar nerve (5). Physiological sodium chloride was used as a cooling solution. For all osteotomies mode boosted burst c and pump 5 were used.

Results: Piezoelectric osteotomy permitted micrometric selective cutting and a clear surgical site due to the cavitation effect created by the cooling solution and the oscillating tip. No excessive bleeding was encountered. The risk of accidental soft tissue harm, such as perforating the sinus membrane or damaging adjacent nerves, was definitely lower than in the case of a conventional bur. No serious complications were encountered in the postoperative wound healing process after 2, 14, 30 and 90 days. However, the surgical procedures were time-consuming.

Conclusion: Piezosurgery is an advantageous osteotomy technique for delicate structures in the oral and maxillofacial region. With respect to osteotomies of thin and fragile bones, the application of ultrasound is superior to other mechanical instruments because of the extremely precise and virtually arbitrary cut geometries, easy handling, efficient bone ablation and minimal accidental damage to adjacent soft tissue structures.

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ENGLISH Br J Oral Maxillofac Surg. 2008; 46(4):265-269

■ Piezoelectric surgery: Twenty years of use.
 The use of ultrasonic vibrations for the cutting of bone was first introduced two decades ago. Piezoelectric surgery is a minimally invasive technique that lessens the risk of damage to surrounding soft tissues and important structures such as nerves, vessels, and mucosa. It also reduces damage to osteocytes and permits good survival of bony cells during harvesting of bone. Piezoelectric surgery was first used by oral and maxillofacial surgeons for osteotomies, but recently some specific applications in neurosurgery and orthopaedics have been proposed. We review the different applications of piezoelectric surgery.

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ENGLISH J Oral Maxillofac Surg. 2005; 63(9):1283-1287

■ Intraoral Piezosurgery: Preliminary Results of a New Technique.
 The piezosurgery instrument, developed in 1988, uses a modulated ultrasonic frequency that permits highly precise and safe cutting of hard tissue. Nerves, vessels, and soft tissue are not injured by the microvibrations (60 to 200 mm/sec), which are optimally adjusted to target only mineralized tissue. The selective and thermally harmless nature of the piezosurgery instrument results in a low bleeding tendency. In addition, the instrument can be used in operations requiring either local or general anaesthesia. The precise nature of the instrument allows exact, clean, and smooth cut geometries during surgery. The difference in time requirement for surgical procedures using the piezosurgery instrument in comparison with the conventional drill is negligible. Postoperatively, excellent wound healing, with no nerve and soft tissue injuries, is observed. It is apparent that the range of application of piezosurgery is not limited to minor operations. Because of its highly selective and accurate nature, with its cutting effect exclusively targeting hard tissue, its use may be extended to more complex oral surgery cases, as well as to other interdisciplinary problems.

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